

Electrical Conductivity of Aqueous Solutions

PRE-LAB ASSIGNMENT:

Reading: Chapter 4.1-4.3 in Brown, LeMay, Bursten & Murphy.

1. Using Table 1 in this handout, determine which solution has a higher conductivity, 0.1 M HCl or 0.1 M NaCl. Explain what creates the higher conductivity.
2. Write the balanced complete ionic equation for the titration of aq. NaOH with aq. HCl.

INTRODUCTION: In this lab you will explore the nature of aqueous solutions by investigating the relationship between conductivity and strong and weak electrolytes. To do this, you will add increasing amounts of either acid or base to several electrolyte solutions. After each addition you will measure the conductivity of the solution. From the conductivity data, you will work to deduce the nature of the chemical reaction that occurred in the experiment. This experiment is another "discovery" exercise to get you thinking about how the experimental data you have collected relates to the concepts you are learning in class. You will work in small groups to decide how best to work up the data. Plotting your data will allow you to visually follow the chemical reactions that are occurring and deduce what ions that are present.

What is conductivity?

Conductivity is a measure of the concentration of ions in solution. By completing the circuit shown in Figure 1, we can measure the conductivity of the solution in the beaker. The conductivity is proportional to the current that flows between the electrodes. For current to flow, ions must be present in solution to carry the charge from one electrode to another. Increasing the number of ions in solution will increase the amount of charge that can be carried between electrodes and will increase the conductivity. The units microSiemens/cm (uS/cm) and milliSiemens/cm (mS/cm) are most commonly used to describe the conductivity of aqueous solutions.¹




Figure 1. Schematic of a simple conductivity measurement system.

¹ The Siemen was formerly called mho (pronounced "mo"), which was derived as a unit of conductivity by reversing the letters in "ohm," the unit of resistance.


Another factor in conductivity measurements is that not all ions carry charge (conduct electricity) equally. H^+ and OH^- are unique and move through solution very rapidly and are very good charge carriers. Ions such as NH_4^+ and Cl^- move through solution at a slower rate and therefore do not conduct electricity as well. Uncharged species in solution do not carry any charge. Table 1 is a table of molar conductivity for the ions in this exercise. The molar conductivity is the conductivity of a solution for the ion containing one mole of charge per liter. Note that the molar conductivity of H^+ ions is 5-7 times the conductivity of other small cations. The molar conductivity of OH^- is 3-5 times the conductivity of other small anions. To calculate the conductivity of a solution you simply multiply the concentration of each ion in solution by its molar conductivity and charge then add these values for all ions in solution.

Table 1. Molar Conductivity of Selected Ions

Ion	Molar Conductivity ($\text{S L mol}^{-1} \text{cm}^{-1}$)
<u>Cations:</u> H^+	0.34982
Na^+	0.05011
NH_4^+	0.0735
K^+	0.0735
<u>Anions:</u> OH^-	0.1986
Cl^-	0.07635
CH_3COO^- *	0.0409
Br^-	0.0781

* CH_3COO^- is the acetate ion, and it is often written as $\text{C}_2\text{H}_3\text{O}_2^-$.

How can we use conductivity to study solutions? Substances that completely dissociate into ions (strong electrolytes) produce solutions with high conductivity. Substances that partially dissociate (weak electrolytes) will produce solutions with low conductivity. Substances that dissolve but do not dissociate produce solutions that have the very low conductivity of pure water. How can we use conductivity to study chemical reactions? Since conductivity is a function of both the concentration and the composition of the solution being measured, we can use conductivity to follow chemical reactions. Consider an acid base titration of a solution of KOH with HBr. The reaction is represented in the following:



During the course of the titration OH^- is consumed, Figure 2. The equivalence point is reached when the number of moles of H^+ added is equal to the number of moles of OH^- in the starting solution. After the equivalence point, the moles of H^+ increase, because they are no longer consumed by reaction with OH^- . For the spectator ions, the number of moles of K^+ remains constant and the number of moles of Br^- increases as the titrant is added. How will these ion concentrations be reflected in the conductivity of the solution?

Initially the solution conductivity is high, because OH^- has a large molar conductivity. However, as OH^- is neutralized the conductivity falls. After all the OH^- is neutralized the H^+ concentration increases and the conductivity again rises rapidly. However, remember that the conductivity includes contributions from all ions in solution, not just H^+ and OH^- .




Figure 2. Titration of 150 mL of 0.00667 M KOH with 0.02 M HBr. Note the expanded scale for Br⁻.

PROCEDURE: For this experiment you will work in pairs and perform conductivity titrations on three different solutions. Every pair will perform a titration of NaOH with HCl. Half of the class will also titrate NH₃ with HCl and NH₄Cl with NaOH. The other half of the class will titrate acetic acid (CH₃COOH) with NaOH and sodium acetate (NaC₂H₃O₂) with HCl. At the end of class you will work in groups of four to analyze your data so that each group has all five of the titrations plotted. You will use this conductivity data to predict the ionic composition of your test solutions before, during, and after the titrations.

Measuring conductivity: You will be using a pH meter in conductivity mode (model: AR20) with a conductivity electrode to analyze your solutions. The meter will be setup for conductivity measurement, so that you will not need to press any buttons on the meter. The electrode will be standardized prior to lab time.

When using the electrode, the level of aqueous solution should be in, but not over, the “circular white zone.” The “white zone” allows the solution to be accurately measured (the solution must cover the “hole” in the electrode body so that sample can reach the electrode); while keeping the electronics dry. The electrode can be easily raised and lowered to keep the solution in this zone. Read the conductivity on the display when the meter reads “stable.” Remember to record the unit of conductivity that the meter is displaying, as the meter will automatically change from uS/cm to mS/cm. When not in use, the electrode is stored in a small beaker of distilled water.

Rinse the electrode with distilled water and dry it gently with a Kim Wipe before each use. This will clean the electrode and minimize contamination between solutions. Use the “waste beaker” to collect the rinse water. This and all other solutions may be poured down the drain (but not the stir bar!).

You and your partner will be performing three titrations together. Table 2 outlines the starting solution and titrant for each group. Notice that there are two different concentrations of NaOH present in the lab; be careful to use the correct one for your titrations.

Table 2. Solutions to be titrated:

<u>Group</u>	<u>Starting Solution</u>	<u>Titrant</u>
Group 1	150 mL of 0.00667 M NaOH	0.02 M HCl
	150 mL of 0.00667 M NaC ₂ H ₃ O ₂	0.02 M HCl
	150 mL of 0.00667 M HC ₂ H ₃ O ₂	0.02 M NaOH
Group 2	150 mL of 0.00667 M NaOH	0.02 M HCl
	150 mL of 0.00667 M NH ₃	0.02 M HCl
	150 mL of 0.00667 M NH ₄ Cl	0.02 M NaOH

Titration procedures for this lab: Place 150 mL of the starting solution in a clean 250-mL beaker. The graduations on your beaker are of sufficient accuracy to measure the solution volume for this experiment. Add a stir bar to the beaker and stir the solution slowly using a magnetic stirrer. Do not allow the stir bar to hit the conductivity electrode. Measure the conductivity of the starting solution.

Titrate the starting solution by adding 10 mL increments of the titrant to the solution (use a graduated cylinder to measure the titrant volume). Adjust the probe height to ensure that the liquid surface is within the probe’s white zone. Allow the solution to mix for a few seconds after each addition before measuring the conductivity. Record your results in your lab book. Add a total of 90 mL of titrant to the beaker. Rinse your beaker and electrode well with distilled water before the next titration.

DATA ANALYSIS: To analyze data for all 5 titrations, you and your partner will need to share data with another student pair. If you collected Group 1 data, you’ll need to get together with some other pair that collected Group 2 data. Make a suitable plot of the conductivity for each titration. Each group of 2 students will need to make its own graph of the five titrations.

What should be in your lab notebook?

Keep track of all measurements and *appropriate units*. Record all data and observations as they occur. Record the names of all the people who contributed data for your group graph.

What should be attached to your discussion now but put into your lab notebook later?

All five final graphs (one per titration) or one graph of all 5 titrations.

The Excel spreadsheet data used to make your graph of all 5 titrations.

A balanced complete ionic equation for each of the 5 titrations (include states of matter & charges).

The completed Supplemental Report Sheet.

What should be in your discussion?

Start off by stating what kind of electrolyte each of the five starting solutions is (a strong or weak electrolyte) and specifically why.

Explain why the conductivity of the solutions during the titrations never reaches zero (this explanation need only be a part of one of your titration analyses).

In your discussion, explain the shape of **EACH** plotted titration graph in its own paragraph using the following parameters as your guideline.

- Before titration, where is the starting conductivity and why?
- Does the conductivity increase or decrease upon initial addition of titrant?
- Chemically, why is the slope increasing or decreasing?
- On the graph there is an inflection point, at what point does that occur?
- Chemically what has changed at this particular point?
- Does the conductivity increase or decrease after the inflection point?
- Chemically, why is the slope increasing or decreasing past the inflection point?

Your discussion should consist of 5 short paragraphs – one for each titration. The one page limit does not need to be adhered to this time, but you should still try to be as brief as possible.

Supplemental Report Sheet

(In each box, fill in all the chemical species present)

Starting Solution/ Titrant	0 mL titrant added	50 mL titrant added	80 mL titrant added
NaOH / HCl			
NaC₂H₃O₂ / HCl			
HC₂H₃O₂ / NaOH			
NH₃ / HCl			
NH₄Cl / NaOH			