

Lab 7
2-color flow cytometry

Prior to lab, cells will be harvested from a mouse spleen. Red blood cells will be lysed, and the leukocytes resuspended in staining media (SM), ready for you to stain with fluorescent-labeled antibodies.

Stains available

FL	PE
anti mouse CD4-FL	anti-mouse CD3-PE
anti-mouse B220-FL	anti mouse CD8-PE
anti-mouse IgM-FL	anti-mouse CD69-PE
anti-mouse CD69-FL	

Prior to lab:

Think of a simple question you'd like to answer using two of the reagents above. (Don't worry, this can be a VERY simple question)

Question:

Two stains to be used (remember, one must be FL-labeled, and the other PE-labeled):

2-Color Flow Cytometry

Three important things to keep in mind while staining cells for flow cytometric analysis:

- a. keep the cells on ice.
- b. keep stained cells in the dark (under foil) when possible so fluorescent molecules on the antibodies don't fade. The foil can go over the top of the whole ice bucket.
- c. be sure that the final prep is a single-cell suspension, or you risk clogging the intake tube on the flow cytometer. We'll filter the cells through a 100 μm mesh today to ensure this.

Procedure

Mouse leukocyte suspensions were prepared from spleens prior to class.

1. place 100 μl staining media (SM) in a labeled FACS tube
2. add 2 μl of each of the two appropriate fluorescent antibodies to this SM, making sure that the 2 μl of stain is actually delivered into the SM (not stuck on the interior wall of the tube)
 - * do not use the same pipet tip for the 2 Abs – use a NEW tip each time
3. add 100 μl spleen suspension to the stain mixture; tap tube or pipet to mix well
4. place tube on ice, covered with foil, for 20 minutes.
5. add 200 μl SM, then centrifuge 3 minutes at 300 xg.
6. remove supernatant by pipetting, tap tube to loosen pellet (pellets may appear as a solid white circle or just a filmy ring, depending on the cell numbers)
7. resuspend pellet in 400 μl SM, spin in centrifuge; repeat wash (steps 6, 7) one more time
8. resuspend in 1 ml SM
9. filter cells through 100 μm Nitex mesh square into a new FACS tube, using a glass pasteur pipette. Keep covered on ice until your turn to use the flow cytometer.
10. analyze on flow cytometer

* we will hold a lottery during the staining procedure to determine the order in which the samples will be analyzed